

Exploring and Distilling Posterior and Prior Knowledge for Radiology Report Generation

Fenglin Liu¹, Xian Wu², Shen Ge², Wei Fan², Yuexian Zou^{1,3}

¹ ADSPLAB, School of ECE, Peking University

² Tencent Medical AI Lab, Beijing, China ³ Peng Cheng Laboratory, Shenzhen, China



Contents

- Introduction
 - Radiology Report Generation
 - Motivations: Visual Data Deviation
 - Motivations : Textual Data Deviation
- Approach: Posterior-and-Prior Knowledge Exploring-and-Distilling
 - Posterior Knowledge Explorer
 - Prior Knowledge Explorer
 - Multi-domain Knowledge Distiller
- Experiments
 - Quantitative Results
 - Qualitative Results
- Conclusions




1. Introduction



Radiology Report Generation

- Dataset: (V, S) , where V and $S = \{s_1, s_2, \dots, s_T\}$ represent the **input radiology image** and the **target report**, respectively.
- Encoder-Decoder Framework: In the encoding stage, the **global image features** are extracted by CNN from the entire image; In the decoding stage, the **whole report is generated using HRNN**.
- Training Objective: The widely-used training objective is to **minimize the cross entropy loss**.




Visual Enc. : $\mathcal{V} \rightarrow \hat{\mathcal{V}}$; Target Dec. : $\hat{\mathcal{V}} \rightarrow \mathcal{S}$.

$$L_{CE}(\theta) = - \sum_{t=1}^T \log \left(p_{\theta} \left(s_t^* | s_{1:t-1}^*; \mathcal{V} \right) \right)$$



Motivations: Visual Data Deviation



Medical Report:

Lungs are clear. No pleural effusions or pneumothoraces. Heart and mediastinum of normal size and contour. scoliosis.

- The normal images **dominate** the dataset over the abnormal ones [1], especially for the rare diseases.
- For each abnormal image, the appearance of abnormal regions (**red** bounding box) only **occupy a small part** of the entire image.
- As a result, this **unbalanced** visual distribution would **distract** the model from accurately capturing the features of rare and diverse abnormal regions




Motivations: Textual Data Deviation



Medical Report:

Lungs are clear. No pleural effusions or pneumothoraces. Heart and mediastinum of normal size and contour. scoliosis.



Medical Report:

The heart size is enlarged. The aorta is tortuous. The pulmonary vasculature appears normal. Lungs are otherwise clear bilaterally. No pleural effusions or pneumothorax. No bony abnormalities.

- In a report, radiologists tend to describe all the items in an image, making the description of abnormal regions (red colored text) only occupy a small part of the entire report.
- Besides, there are many similar sentences (blue colored text) used in each report to describe the normal regions.
- With this unbalanced textual distribution, training with such dataset makes the generation of normal sentences dominant, disabling the model to describe specific crucial abnormalities.



2. Approach: Posterior-and-Prior Knowledge Exploring-and-Distilling




Overview: PPKED

- The Posterior-and-Prior Knowledge Exploring-and-Distilling (PPKED) imitates the **radiologists' working patterns**. Given a medical image, radiologists **will**:
 - 1. **examine the abnormal regions and assign the disease topic tags to the abnormal regions**;
 - 2. then accurately write a corresponding report based on years of **prior medical knowledge** and **prior working experience** accumulations.



Overview: PPKED

To model above working patterns, the PPKED introduces Posterior Knowledge Explorer (**PoKE**), Prior Knowledge Explorer (**PrKE**) and Multi-domain Knowledge Distiller (**MKD**).




Our approach based on the Multi-Head Attention (**MHA**) and Feed-Forward Network (**FFN**) from Transformer [1].


[1] Attention is all you need. In NIPS, 2017.



Posterior Knowledge Explorer (PoKE)



Posterior Knowledge Explorer (PoKE)



Given the input image I and disease topics tags T :


$$\hat{T} = \text{FFN}(\text{MHA}(I, T)); \hat{I} = \text{FFN}(\text{MHA}(\hat{T}, I))$$

$$I' = \text{LayerNorm}(\hat{I} + \hat{T})$$

i.e., the I are first used to **find** the most relevant topics and filter out the irrelevant topics. Then the attended topics \hat{T} are further used to **mine** topic related image features \hat{I} .




Prior Knowledge Explorer (PrKE)



The PrKE is designed to **alleviate textual data deviation** by encoding the prior knowledge.



Prior Knowledge Explorer (PrKE)



The prior knowledge includes the **prior radiology reports** W_{Pr} (i.e., prior working experience) pre-retrieved from the training corpus and the **prior medical knowledge graph** G_{Pr} (i.e., prior medical knowledge), which models the domain-specific prior knowledge structure and is pre-constructed from the training corpus:


$$W'_{Pr} = \text{FFN}(\text{MHA}(I', W_{Pr}))$$

$$G'_{Pr} = \text{FFN}(\text{MHA}(I', G_{Pr}))$$

By processing I' through these two equations, we can acquire W'_{Pr} and G'_{Pr} which represent the prior knowledge **relating** to the **abnormal regions** I' of the input image.




Knowledge Distiller (MKD)



Finally, the MKD is designed to **distill** the useful knowledge to generate proper reports.



Knowledge Distiller (MKD)



Given the embedding of current input word x_t :

$$h_t = \text{MHA}(x_t, x_{1:t})$$

$$h'_t = \text{ADA}(h_t, I', G'_{\text{Pr}}, W'_{\text{Pr}})$$

$$\text{ADA}(h_t, I', G'_{\text{Pr}}, W'_{\text{Pr}}) = \text{MHA}(h_t, I' + \lambda_1 G'_{\text{Pr}} + \lambda_2 W'_{\text{Pr}})$$

$$\lambda_1, \lambda_2 = \sigma(h_t W_h \oplus (I' W_I + G'_{\text{Pr}} W_G + W'_{\text{Pr}} W_W))$$

$$y_t \sim p_t = \text{softmax}(\text{FFN}(h'_t) W_p + b_p)$$

where x_t denotes the embedding of current input word; y_t denotes the current target word; σ and \oplus denote the sigmoid function and the matrix-vector addition, respectively; ADA denotes the Adaptive Distilling Attention; The λ_1 and λ_2 weight the **importance** of G'_{Pr} and W'_{Pr} for each target word, respectively.



3. Experiments



Quantitative Results

Dataset	Methods	Year	BLEU-1	BLEU-2	BLEU-3	BLEU-4	METEOR	ROUGE-L	CIDEr
MIMIC-CXR	CNN-RNN	2015	0.299	0.184	0.121	0.084	0.124	0.263	-
	AdaAtt	2017	0.299	0.185	0.124	0.088	0.118	0.266	-
	Att2in	2017	0.325	0.203	0.136	0.096	0.134	0.276	-
	Up-Down	2018	0.317	0.195	0.130	0.092	0.128	0.267	-
	Transformer	2020	0.314	0.192	0.127	0.090	0.125	0.265	-
	R2Gen	2020	0.353	0.218	0.145	0.103	0.142	0.277	-
	PPKED	Ours	0.360	0.224	0.149	0.106	0.149	0.284	0.237
IU-Xray	HRNN	2017	0.439	0.281	0.190	0.133	-	0.342	0.261
	CoAtt	2018	0.455	0.288	0.205	0.154	-	0.369	0.277
	HRGR-Agent	2018	0.438	0.298	0.208	0.151	-	0.322	0.343
	CMAS-RL	2019	0.464	0.301	0.210	0.154	-	0.362	0.275
	Transformer	2020	0.396	0.254	0.179	0.135	0.164	0.342	-
	R2Gen	2020	0.470	0.304	0.219	0.165	0.187	0.371	-
	PPKED	Ours	0.483	0.315	0.224	0.168	0.190	0.376	0.351

Table 1. Results of the PPKED and other methods on MIMIC-CXR [1] and IU-Xray [2] datasets.

[1] MIMIC-CXR: A large publicly available database of labeled chest radiographs. arXiv preprint arXiv:1901.07042, 2019.

[2] Preparing a collection of radiology examinations for distribution and retrieval. Journal of the American Medical Informatics Association, 23(2):304–310, 2016.



Qualitative Results

Figure 1. Two examples of ground truth reports and reports generated by HRNN [1] and our method. The **Red** colored text indicates the abnormalities. The **HRNN fails** to depict some rare but important abnormalities and generates some error sentences (**Blue** colored text) and repeated sentences (**Underlined** text). **Our PPKED has higher rate of accurately describing the rare and diverse abnormalities.**



Ground Truth:

Lungs are clear. No pleural effusions or pneumothoraces. Heart and mediastinum of normal size and contour. ¹**scoliosis**.

HRNN:

Heart size is normal. There is a moderate right sided pneumothorax with tip in the right atrium. There is a moderate right sided pneumothorax with large pleural effusion. No pneumothorax masses. No pneumothorax masses. No acute bony abnormalities.

Ours:

¹**There is a scoliosis.** No acute cardiopulmonary abnormality. There is no pleural effusion. No evidence of pneumothorax. The lungs are clear. There is no focal airspace consolidation.



Ground Truth:

¹**The heart size is enlarged.** ²**The aorta is tortuous.** The pulmonary vasculature appears normal. Lungs are otherwise clear bilaterally. No pleural effusions or pneumothorax. No bony abnormalities.

HRNN:

¹**Cardiomegaly with pulmonary vascular congestion and interstitial edema.** There is a moderate right sided pneumothorax with large pleural effusion. No bony abnormalities. There is no pneumothorax. There is no pneumothorax.

Ours:

¹**Heart size is enlarged.** ²**Tortuosity of the aorta.** No pleural effusion. There is no focal airspace consolidation. There is no pneumothorax. No bony abnormalities.



4. Conclusions



Conclusions

- In this work, we present an effective approach of **exploring and distilling posterior and prior knowledge** for radiology report generation.
- Our approach **imitates** the working patterns of radiologists to **alleviate** the data bias problem.
- The experiments demonstrate the **effectiveness** of our method.
- Our approach not only generates **meaningful** and **robust** radiology reports supported with **accurate** abnormal descriptions and regions, but also **outperforms** previous state-of-the-art models on the two public datasets.





Thank you for your attention!

If you have any questions about our paper, you can send an email to fenglinliu98@pku.edu.cn

